

# The SciLog<sup>®</sup> LabTec<sup>™</sup> is an automated laboratory scale dispensing system by volume, weight or weight ratio.

The LabTec<sup>™</sup> automates, optimizes and documents repetitive liquid dispensing with or without in-line filter sterilization. The LabTec<sup>™</sup> models are ideally suited for dispensing sterile solutions in media kitchens and microbiology laboratories. Also widely used in small production runs to fill your final product into vials, bottles, bags, and containers.

The automatic documentation and alarm / pump stop settings allow the user to focus on other tasks while the system is running. Remote control and programmable end points ensure the system starts and stops operating when a given dispense is complete or an interlock condition occurs. When sold with SciDoc software or a printer, documentation capabilities include 7 real-time parameters.

#### **Features and Benefits**

- Rapid, high prevision dispensing by volume, weight or weight ratio
- Self priming
- Reversible flow
- Remote activation
- Accuracy +/- 0.5% by volume

- Accuracy +/- 0.3% by weight
- Pressure alarm alerts when in-line filter needs to be replaced
- PC or printer documentation of dispensing

## SciLog<sup>®</sup> LabTec<sup>™</sup>

- intelligent bioprocessing system
- dispensing system



### **Performance Characteristics**

Table 1 - Automated Weighing and Dilution of Food Samples						
Trial	Sample Weight (g)	Theoretical Diluent Weight (g)	Actual Diluent Weight (g)	Errors (g)		
1	10.0	90.0	90.0	0.0		
2	10.0	90.0	89.9	-0.1		
3	10.0	90.0	90.1	+0.1		
4	10.0	90.0	89.8	-0.2		
5	10.0	90.0	90.1	+0.1		
6	10.0	90.0	90.0	+0.1		
7	10.0	90.0	89.9	-0.1		
8	10.0	90.0	90.2	+0.2		
9	10.0	90.0	90.1	+0.1		
10	10.0	90.0	89.8	-0.2		

Table 2 - High Speed Volumetric Dispensing						
Dispensed Volume (ml)	Tubing Size	Pump Speed (%)	Slow Factor (ml)	Typical Precision (%)	Dispensing Time / Aliquot (sec)	
5.0	#15	100%	2.50	1.7%	2.1	
10.0	#15	100%	2.50	0.5%	2.0	
25.0	#15	100%	2.50	0.20%	2.7	
25.0	#24	100%	4.75	0.15%	2.9	
50.0	#24	100%	4.75	0.50%	3.9	
100.0	#24	100%	4.75	0.35%	6.3	
100.0	#35	100%	5.00	<0.5%	4.5	
200.0	#35	100%	5.00	<0.5%	6.0	
300.0	#35	#35 100%		5.00 <0.5%		




Weight Entered (g)	Dispensed Weight (g)	RSD (%)	Time per Aliquot (sec)
200.00	199.95	0.03%	16
150.00	150.01	0.11%	14
100.00	100.01	0.11%	13
50.00	49.96	0.18%	12
25.00	25.05	0.25%	10

### SciLog<sup>®</sup> LabTec<sup>™</sup>


## Applications

#### Weighing, Dilution and Sterilization




The LabTec™ significantly increases the productivity of any microbiology laboratory. The dispensing system automates weighing and diluting of food samples, in preparation for bacteriological analysis. Only an approximate sample size is needed. The LabTec™ dispenser calculates and rapidly dispenses the media required to achieve a user-selected diluent factor. See Table 1. Costly and timeconsuming autoclaving of media is avoided with in-line filter sterilization of diluents using a filter capsule. A disposable pressure sensor monitors filter backpressure and provides an alarm signal when a failing sterilizing filter (high backpressure) needs to be replaced. An optional printer or PC hook-up automatically documents all sample and diluent weights. The LabTec™ dispenser improves the speed, precision and accuracy of the entire sample preparation process by automating the weighing, dilution and documentation process.

#### Precision Batch Dispensing



The LabTec<sup>™</sup> Smart Dispensing System is capable of high speed, high precision batch dispensing. The LabTec<sup>™</sup> comes with a 1082 Tandem peristattic pump head. Up to 10 different dispensing volumes can be stored and easily retrieved for quick batch volume dispensing. An optional sterilizing filter and disposable pressure sensor provide for in-line buffer sterilization and reliable detection of filter plug-up conditions. Typical performance results including dispensing times are summarized in Table 2. The LabTec<sup>™</sup> is easy to calibrate. A stored calibration curve is provided for each pump head/pump tubing combination. From a stored menu, you select the pump tubing you have installed in the LabTec<sup>™</sup>. For high accuracy dispensing applications, you may want to use the LabTec<sup>™</sup>'s single-point re-cal feature.

#### Fast, Accurate Dispensing and Sterilization



In this configuration, the LabTec™ is connected to an electronic scale for high accuracy filling applications. A sterilizing filter and a disposable pressure transducer provide a safe and effective in-line filter sterilization capability. The LabTec™ continuously monitors the filter backpressure and alarms when a user-defined pressure level has been exceeded, which indicates a filter plug-up condition. See Table 3. The high dispensing accuracy is achieved by reducing the pump rate as final target weight is being approached. The slow-down avoids overshooting the target weight. After slow-down, the LabTec™ pump stops briefly when 99% of the final target weight has been dispensed. The electronic scale is allowed to come to a steady-state readout and the LabTec™ slowly starts up again to dispense the remaining solution.

### **Specifications**

	Description
Dimension / Weight	Width: 5.75" (146 mm) x height: 8.5" (2126 mm) x depth: 11" (279 mm): 14 lbs (6.4 kg)
Enclosure & Rating	16 Ga, aluminium baked epoxy blue 4-40dC, 0-100% humidity, IP20
Pressure Sensors	Accommodates one disposable pressure sensor. The calibrated pressure range is 0 - 60 psi. Any point within this range can be recalibrated using an external pressure reference source.
Power	115 / 220-240 VAC, 60 / 50 Hz, 75 Watts, double fused: T1AL 250V (CE: IR35A 250VAC)
Motor / Encoder 8, 160, 600, 3400 RPM, 30 VDC, 3.8A, 100 ppr	
I/O Ports	Male DB9 scale connections (RS-232), female DB9 printer or PC connection (RS-232), external IO DB37 connector, 1 TTL input, 4 TTL output, 3 4-20mA
Operational Mode	Dispense by volume, weight or weight ratio. Store up to 10 programs per mode. Pump re-calibration.

### **Options and Accessories**

#### Pump Heads:

- Tandem Peristaltic
- 1081 flow rate (ml/min): 0.03 1515
- 1082 flow rate (ml/min): 0.5 2258
  (Pressure: 25 psi continuous 45 psi max)
- FMI RH piston: 0.002-320 (600, 3400 RPM) (Pressure: 100 psi max)
- Micropump MAG: 0.54-3488 (3400 RPM) (Pressure: 40 - 70 psi max[model dependant])
- Masterflex peristaltic:
  0.03- 2900 (8, 160. 600 RPM)
  (Pressure: 25 psi continuous 45 psi max)

<b>J</b>	Ordering	Information
----------	----------	-------------

0	- LABT-1									
	Code   Motor			Pump	Head					
	0	8 RPM 160 RPM	Code   SciLog	Code   Masterflex	Code   Micro Pum	р	Code   FMI RH			
	1 3 6	3400 RPM 600 RPM	81 1081 Pump 82 1082 Pump	21 Thin Wall 22 Thick Wall 62 Hi Performance L/S	31 MAG 120 32 MAG 201 33 MAG 1840	35 MAG 200 36 MAG 201 37 MAG 187	41 OCKC 42 1CKC			
					34 MAG 040					

Example: 200-LABT-1181 - SciLog® LabTec™ - 160 RPM motor and 1081 head

Parker domnick hunter has a continuous policy of product development and although the Company reserves the right to change specifications, it attempts to keep customers informed of any alterations. This publication is for general information only and customers are requested to contact our Process Filtration Sales Department for detailed information and advice on a products suitability for specific applications. All products are sold subject to the company's Standard conditions of sale.